INTRODUCTION

Portable massage gun with 6 different speeds and 4 different massage heads.
Please read the manual carefully to understand how to use the Massager.

PARTS INFORMATION

1. Ball head (Designed for massaging large muscle groups)
2. Fork head (Designed for massaging the spine)
3. Flat head (Designed for massaging all parts of the body)
4. Bullet head (Designed for massaging joints, deep tissue and trigger points)
5. Main body
6. Level indicator
7. Level adjustment switch
8. Charging port

LEVEL ADJUSTMENT

Repeatedly press the level switch to increase the level of the massage gun. After level 6 the massage gun will turn into the pause mode.

Level 1	-	2.100 hits / minute
Level 2	-	2.400 hits / minute
Level 3	-	2.700 hits / minute
Level 4	-	3.000 hits / minute
Level 5	-	3.300 hits / minute
Level 6	-	3.600 hits / minute

HOW TO USE?

- Choose the massage head depending on your preference and press it in.
- Press the level adjustment button to start the massager. Choose the level of your preference and press the massager on the part of the body you want to massage.
- After your massage, turn off the massager and take out the massage head.

KEEP IN MIND

- Use the massage gun not longer then 15 minutes in one go to protect it for overheating. The massage gun will automatically turn off if the temperature gets too high.
- The massage gun can not be used when charging. Make sure to fully charge the massage gun each time.
- The massage gun can only be used on a dry surface.

WARNING

In order to reduce the potential risk of shock, fire or injury, the product must be used via the following instructions:

- The massage gun can only be used by adults.
- Do not massage one area for more then 1 minute. Only use the massage gun on a dry surface. Massaging same position for along te may casue abrasion of the body.
- Always stop when feeling pain or discomfort. The high frequency levels can only be used on soft tissue parts.
- Choose the appropriate massage head according to the body part described in the manual.
- Keep your fingers or hair away from the expensaion rod and the main body to avoid potential injury.
- Do not put anything into the heat emission holes on the backside.
- Do not immerse the massage gun in water and avoid water in to the device.
- Do not discard or misuse this device.
- Only use the provided charger to charge the device.
- Please check the evice carefully before use.
- Do not make any modification to the device privately.
- Do not run or charge the device unmanned.

CLEANING THE DEVICE

- Make sure the device is turned off before you start cleaning the device.
- Wipe the surface of the device with a slightly damp towel, then dry it after with soft cloth.
- Clean the device after each usage carefully.

TECHNICAL SPECIFICATIONS

Product weight:	935 gram
Material:	ABS + TPR + PP + Zinc Alloy + EVA
Input:	110-240 volt
Battery capacity:	2000 mAh
Rotate speed:	2.100 ~ 3.600 hits / minute
Power:	20 watt
Color:	Silver / Space Grey
Frequencies:	6-level variable frequency vibration

© 2021 The Investment Company Europe B.V. Hyundai is a trademark of Hyundai Corporation Holdings, Korea. Licensed to The Investment Company Europe B.V. All rights reserved.

Licensed by Hyundai Corporation Holdings, Korea. Imported/Distributed by The Investment Company Europe B.V. exclusively developed for sales in the Netherlands, Belgium, and Luxembourg.

PATENT PENDING

The artwork and design of this package are protected by European copyright law and may not be reproduced, distributed, displayed, published or used for any purpose without prior written permission. It is not permitted to alter, or reproduce any trademarks or copyright notice from this package.

KEEP THIS INSTRUCTION MANUAL FOR REFERENCE AS IT CONTAINS IMPORTANT INFORMATION.

FOR QUESTIONS OR SUPPORT PLEASE CONTACT US AT:

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